

# Thermal Vacuum Chamber

[csum-contact@umontpellier.fr](mailto:csum-contact@umontpellier.fr)



Radiative Temperature Interface

GN2 Screen

Maximum

+150°C

Minimum

-170°C

Conductive Temperature Interface

Table

Maximum

+70°C

Minimum

-70°C

No-load ultimate vacuum

$P < 1.10^{-5}$  mbar

Maximum temperature transition speed

$\pm 2^\circ\text{C} / \text{minute}$

Recommended Specimen size

50 x 50 x 50 cm

Maximum Specimen size

70 x 80 x 60 cm

Maximum Specimen mass

50 Kg

T-QCM



ISO 8 Clean Room

